REF V25

Pipette 25ml

25ml

10ml

25ml

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Operation of the instrument

The 25 ml volume is completely dispensed out of the burette when you start with the red line exactly at the 0 graduation and titrate until the red line is exactly at the 25 ml mark.

Hold the instrument in the palm of your hand with the piston rod extending from the top under your thumb. Roll the piston rod with your thumb (or pull the small ball at the end of the piston rod):

- upward to draw in the liquid,
 - or downward to dispense it.
- Fill the burette smoothly and consistently.

Before filling, wet the seal with a small amount of the liquid being used by drawing and pushing the liquid back and forth with your thumb in the area of the burette tip.

Begin filling the burette with the conical sealing joint fully inserted into the bottom of the burette tip.

Start titration with the burette tip touching the side of the dosing container.

Important note

It does not matter if there is a small air bubble below the conical seal during filling, as long as the bubble is not dispensed.

Cleaning

For immediate reuse with another chemical product, fill and empty the burette several times with water. Pump the seal rapidly in the area of the burette tip under running water to facilitate cleaning.

Take a small amount of the new chemical liquid you will use and pump up and down to wet the entire internal surface of the burette, then discard this liquid.

Note

For greater precision during titration, all types of pipettes/burettes should be pre-conditioned with the new fluid to be used.

Our instrument makes this step safer, faster, and more practical, while reducing chemical consumption.

- **General use:** Clean the inside of the burette by filling and emptying it with water.
 - **Deep cleaning:** Remove the burette and wash all parts separately.
 - Our instrument is autoclavable at 120°C for 20 minutes. (Preferably remove the piston from the burette beforehand.
-
- It is recommended to keep the instrument parts assembled during cleaning to minimize the need for adjustments later.
 - Keep the sealing area clean, avoiding stagnant chemicals for too long—ensure the burette remains humidified inside.
 - **To increase or decrease the sealing diameter, proceed as follows:**
1. Turn the knurled seal screw clockwise to tighten the conical seal and increase sealing (or counterclockwise to loosen it).

2. Reinsert the seal into the burette and test the result. Repeat if necessary.
 3. Retract the piston rod fully into the handle before reconnecting the burette.
- No leakage should occur as the sealing diameter is adjustable and represents the only sealing point.
 - All parts are replaceable and fully repairable.
 - Turn the knurled seal screw clockwise to tighten the conical seal and increase sealing.
 - This prevents air bubbles from forming at the top of the cone.

Methanol, acetone, concentrated ethyl alcohol, and benzene are incompatible with the plastic of the burette. Similarly, extremely high concentrations of acid can also affect the plastic.